










Face Shield Assembly

XCL993039 - 4/20



Face Shield Spec Sheet



In response to the urgent need for Personal Protective Equipment (PPE), the Ohio Manufacturing Alliance to Fight COVID-19 has rapidly expanded production of a **reusable face shield** designed to protect front line healthcare workers.

Features:

- Easily assembled
- Reusable
- One size fits all; adjustable for fit and comfort
- Protective shield extends above the head and around the face (see photos for reference)
- Developed in partnership with hospital clinicians
- Able to be disinfected with alcohol and using several other common methods

Minimum Order Quantity (MOQ): 100 units, ready to assemble

Part	Description	Sanitation
Head Harness	Polypropylene 7.5" x 5.7" x .79"	Autoclave* Alcohol Bleach
Head Band	Silicone 15.75" x 1.0" x .06"	
Face Shields (three types)	Polycarbonate #3293Q0000 15.5"x 11.25" Thickness: 0.02"	*Limited autoclave use for polycarbonate shield
	PETG #3293N0000 15.5"x 11.25" Thickness: 0.015" - 0.02"	DO NOT Autoclave Alcohol Bleach
	PET #3292P0000 15.5." x 11.25" Thickness: 0.007"	

Other sterilization methods way be compatible but **have not been reviewed by MAGNET** at the time of this release.

For more information, visit www.repurposingproject.com/shield